Continuum[®] Acetabular System

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1. Zimmer Biomet 2017 Continuum Acetabular Hip System global sales data.
2. Melvin, J.S.et al.Early Failures in Total Hip Arthroplasty.J. Arthroplasty 2014; 29:1285-1288
3. Haynes, J. A. et al. Contemporary Surgical Indications and Referral Trends inRevision Total Hip Arthroplasty: A 10-Year Review.J. Arthroplasty 2016; 31:622-625.
4. Sadgohi, P. et al. Revision Surgery after Total Joint Arthroplasty: A Complication Based Analysis Using Worldwide Arthroplasty Registers.J. Arthroplasty 2013; 28:1329-1332.
5. Rajae, S. et al. Increasing Burden of Total Hip Arthroplasty Revisions in Patients Between 45 and 64 Years of Age.JBJS 2018; 100:449-458.
6. Remes, V. et al. Continuum Ceramic Bearing System In Total Hip Arthroplasty - A Multi-Centre Pan-European Study.Paper 2156.16th EFORT Annual Congress.2015年
7. Parvizi, J. et al. Randomized Controlled Trial of Contemporary Cups with Contemporary Bearings: Trabecular Metal and Ceramic. Abstract.2017.
8. Matharu, GS.et al.Trabecular metal acetabular components reduce the risk of revision following primary total hip arthroplasty: A propensity score matched study from the National Joint Registry for England and Wales.Journal of Arthroplasty. 33(2): 447–52, 2018.
9. Xu, et al. Three-Dimensional Host Bone Coverage in Total Hip Arthroplasty for Crowe Types II and III Developmental Dysplasia of the Hip.Journal of Arthroplasty.2016.
10. Taheriazam A., et al. Bilateral Total Hip Arthroplasty in Femoral Head Avascular Necrosis: Functional Outcomes and Complications. International Journal of Medical Research & Health Sciences 2016; 5(6):51-56.
11. Flecher X., et al. Use of porous tantalum components in Paprosky two and three acetabular revision.A minimum five-year follow-up of fifty one hips.International Orthopaedics 2016.
12. Graceffa A., et al. Clinical outcome of design modifications to the CLS Spotorno Stem in total hip replacement.Joints 2016.
13. Yi J., et al. Result of Modular Necks in Primary Total Hip Arthroplasty with a Average Follow-up of Four Years.Hip Pelvis 2016; 28(3):142-147.
14. Macheras, G. et al. Survivorship of a Porous Tantalum Monoblock Acetabular Component in Primary Hip Arthroplasty with a Mean Follow-Up of 18 Years.The Journal of Arthroplasty.32(12): 3680-84, 2017.
15. Sousa, P. et al. Uncemented Highly Porous Tantalum Acetabular Components: What is the Risk of Failure 10 Years After Revision Total Hip.Paper 541.AAOS 2018 Annual Meeting.
16. Tokarski, A., Novack T., Parvizi J. Is tantalum protective against infection in revision total hip arthroplasty? The Bone and Joint Journal.97:45-9, 2015.
17. According to NJR data from 2003 to 2015 where 9,573 Trabecular Metal and 30,452 non-Trabecular Metal cups were used in revision THA and based on hazard ratios adjusted by patient gender, age group, and indications (OA/non-OA). The complete NJR report can be found at: www.zimmerbiomet.com/TM.
18. Teeter, M. et al. Thirteen-year wear rate comparison of highly crosslinked and conventional polyethylene in total hip arthroplasty: long-term follow-up of a prospective randomized controlled trial.Can J Surg.2017 Jun; 60(3): 212–216
19. Bragdon, C. et al. Clinical Multicenter Studies of the Wear Performance of Highly Crosslinked Remelted Polyethylene in THA.Clin Orthop Relat Res (2013) 471:393–402.
20. Currier, et al. Highly Cross-linked UHMWPE Oxidation – an Improvement Over Conventional Gamma-sterilized? Poster 0983, ORS 2014.
21. Mimnaugh, K.et al.100 Million-Cycle Wear Evaluation of Crosslinked Vitamin E Grafted Polyethylene (VE-HXPE) Acetabular Liners.Paper No. 0403.ORS 2016 Annual Meeting
22. Peiserich M, et al. Retention of Mechanical Properties in a Blended Vitamin E Polyethylene After Extreme Oxidative Challenge.Poster 1060, ORS 2013 Meeting.2013.
23. Bragdon, C. et al., Seven to Ten Year Follow-Up of Highly Crosslinked Polyethylene Liners in Total Hip Arthroplasty, Poster No. 2444, 55th Annual Meeting of the Orthopaedic Research Society, Las Vegas, 2009
24. Wannomae KK, et al., In vivo oxidation of retrieved crosslinked ultra-high molecular-weight polyethylene acetabular components with residual free radicals, J Arthroplasty.2006 ; 21(7): 1005-1011
25. Kurtz, et al. Reasons for Revision of First-Generation Highly Cross-Linked Polyethylenes.J Arthroplasty.Vol. 25 No. 6 Suppl. 1 2010
26. Hamilton, William et al. "THA With Delta Ceramic on Ceramic.Results of a Multicenter Investigational Device Exemption Trial."Clin Orthop Relat Res.(2010) 468:358-366.
27. NJR 14TH Annual Report, page 73.
28. Erens, G.A., T.S.Thornhill, and J.N.Katz."Complications of total hip arthroplasty."In.: UpToDate.2017.
29. Y Zhang, et al, Interfacial Frictional Behavior: Cancellous Bone, Cortical Bone, and a Novel Porous Tantalum Biomaterial, Journal of Musculoskeletal Research, 1999 ,3:4, 245-251
30. JD Bobyn, et al, Characterization of a New Porous Tantalum Biomaterial for Reconstructive Orthopaedics, 66th Annual AAOS, 1999.
31. Vanhegan, L. et al. A financial analysis of hip revision arthroplasty. JBJS.94-B, 5, 2012.
32. Zimmer ZRR_WA_2409_11
33. Zimmer ZRR_WA_2401_11, Rev. 2
34. Zimmer ZRR_WA_2413_11
35. Zimmer Technical Memo 1129.98
36. Kreider D, Continuum® and Trilogy® IT Acetabular Systems – Clinical Evaluation Report.Document No. CER 30048.4.0. p 71-72, 4-June-2015.
37. Rajpura, A. et al. "The current state of bearing surfaces in total hip replacement."Bone Joint Journal.2014;96-B:147-156.
38. Zimmer ZRR_WA_3030_15
39. Kuntz M, Validation of a New High Performance Alumina Matrix Composite for use in Total Joint Replacement.Seminars in Arthroplasty 17, 2006: 141-145
40. Hallab, N. J. et al. Differences in the fretting corrosion of metal-metal and ceramic-metal modular junctions of total hip replacements. Journal of Orthopedic Research.2004; 22:250-259.
41. Kurtz, Steven M. et al. Do Ceramic Femoral Heads Reduce Taper Fretting Corrosion in Hip Arthroplasty? A Retrieval Study.Clin Orthop Relat Res.2013.471:3270-3282.
42. Lusardi, A. Zimmer TMT Technical Memo No. 2012-0029. Trabecular Metal Acetabular Components - Bone Ingrowth Claim Documentation.20-Sept-2012.
43. Karageorgiou V, et al., Porosity of 3D biomaterial scaffolds and osteogenesis.Biomaterials.2005; 26: 5474-5491.
44. Testing conducted by independent lab. Report data on file at Zimmer Biomet. Laboratory testing August, 2013. Test results are not necessarily indicative of clinical performance.
45. Latest ODEP ratings can be found at www.odep.org.uk.

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© 2021 Zimmer Biomet